

Customer Experience: Hospital



SITE:

**Turku University Hospital
(TYKS), T-hospital**
Oncology Ward
Turku, Finland



PROBLEM:

At the Oncology Ward there are lots of products being used, which they were unable to properly disinfect. Cooling caps and the cuffs of blood pressure meters, for example, are challenging to disinfect. Interest in the WiSDOM DS surface disinfection box was sparked after hearing that a regional hospital in Salo was using one.

SOLUTION:

For example, the materials used by isolation patients can nowadays be directly taken into the WiSDOM DS for disinfection. The procedure is, that everyone takes their equipment into the device after use. Disinfecting some equipment with other methods can be time consuming, but with WiSDOM DS it is simple and fast – it's finished in just three minutes.

BENEFITS:

"WiSDOM DS is an incredibly easy and handy device – we have already recommended it for our Pulmonary Diseases Outpatient Clinic too.

Our personnel got immediately excited once starting to use the device. It is easy to place objects inside the disinfection box and we use cleaning wipes for removing organic stains.

The surface disinfection box is one more method for reducing the risk of hospital bacteria."

-Paula Halonen, Assistant Head Nurse